
Lime-Lite™ ENHANCED

BIOACTIVE LIGHT CURE CAVITY LINER

XP-LLE-IN-07w
REV: 10/2024

EN

PULPDENT®

LIME-LITE™ ENHANCED

Light Cure Cavity Liner

Calcium, Phosphate, and Fluoride Releasing

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Dentin Shade • Contains no BPA

Lime-Lite Enhanced is a moisture tolerant, light cure, radiopaque cavity liner and base with self-adhesive properties. It is a mineral-enriched urethane dimethacrylate resin that releases calcium, phosphate, and fluoride ions and is formulated for use with adhesives, composites, and conventional restorative materials. It is similar in chemistry to dental composites and chemically bonds to composites and adhesive primers. When increased bond strength to dentin is desired, use of a bonding agent is recommended. Lime-Lite Enhanced contains a patented rubberized-resin component that absorbs shock and stress and resists fracture and chipping. The material sets hard upon light curing, has a dentin shade, and is virtually insoluble in water and oral fluids. Lime-Lite Enhanced contains no BHT, no Bisphenol A, no Bis-GMA, and no BPA derivatives.

PRODUCT DESCRIPTION

Lime-Lite Enhanced is a mineral-enriched urethane dimethacrylate resin cavity liner and base that is moisture friendly, light cure, and radiopaque. It releases calcium, phosphate, and fluoride.

INDICATIONS

Indicated as a light cure cavity liner and base.

CONTRAINDICATIONS

Not indicated for direct placement on the exposed pulp.

INSTRUCTIONS FOR USE

1. Prepare cavity, rinse, and dry, but do not desiccate the tooth.
2. Place a pink 20-gauge dispenser tip securely onto the Lime-Lite Enhanced syringe.
3. Apply 1-1½mm of Lime-Lite Enhanced to the dentin.
4. Light cure for 20 seconds (450 mW/cm²). Do not exceed manufacturer's recommended curing time for the light you are using.
5. If used as a base, apply additional increments of 1-1½mm and light cure each layer. Cap syringe immediately after use.

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Notification des incidents graves

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veuillez consulter le site web [https://www.fda.gov/medical-devices/dental-devices/ multi-ple-use-dental-dispenser-devices](https://www.fda.gov/medical-devices/dental-devices/multi-ple-use-dental-dispenser-devices)

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Ravelli S.p.A.
Hofmeester Dental B.V.
Lifco Dental AB
UNIDENT AB
PULPDENT® Corporation 80 Oakland Street • Watertown, MA 02472 • U.S.A.